

Image Processing Examples

- Example files can be found on the Image Analyst MKII demonstration CD in the folder "Examples", or downloaded from <http://www.imageanalyst.net/tutorials.html>
- The macros used below have been installed into the 'Documents and Settings/username/My Documents/Image Analyst' folder

Spectral Correction, Alignment and advanced FRET

- Click open and select 01a1.nd, 01a2.nd in the "FRET" folder
- In the "File Order" tab of the Multi Dimensional Load dialog make sure, that the order is correctly incrementing
- Switch back to the Open tab, and Click Open
- Select the Visual Macro/New Visual macro menu and load the "5pYC-capain-Fura6F FRET correction.ips" macro
- Run the macro with the blue arrowhead button "M"